

## Minimum technical requirements recommended for the implementation of **NemoStudio Family 2D Software**

<b>Operating System</b>	<b>Minimum:</b> Windows 7 Service Pack 1, Windows 8.x, Windows 10 (professional version); 32 bits. <b>Recommended:</b> Windows 10 Professional 64 bits.
<b>Processor</b>	<b>Minimum:</b> Intel Core i3, i5. <b>Recommended:</b> Intel Core i7 at 2.8Ghz.
<b>RAM Memory</b>	<b>Minimum:</b> 2 Gb. <b>Recommended:</b> 4 Gb.
<b>Graphic Card</b>	<b>Minimum:</b> GPU card with minimum 1 Gb of total VRAM dedicated. <b>Recommended:</b> GPU card with 2 Gb of total dedicated VRAM. (NVIDIA Card is recommended as it is the most tested).
<b>Hard Disc</b>	<b>Between 3 and 5 Gb</b> for the software installation. Additional space may be needed, depending on the volume of cases.
<b>Firewall configuration</b>	Allow access to TCP/IP communication to executables NemoServer.exe and NemoStudio.exe
<b>Network Speed (wired or wireless)</b>	Network local installation. <b>Minimum 1Gb.</b> Internet speed ADSL 50 mb.
<b>Internet Access</b>	Yes, mandatory for licensing control.
<b>Mouse</b>	A mouse with three buttons and a mousewheel is highly recommended.
<b>For Apple Macbook, iMac</b>	<b>Bootcamp (native Windows)</b> for the highest performance. <b>Virtual Pc running Windows</b> in a virtual environment. Parallels is recommended. It is advisable to configure the virtual environment with the recommended Windows requirements aforementioned.