

Minimum technical requirements recommended for the implementation of **NemoStudio Family 3D Software**

Operating System	Minimum: Windows 7 Service Pack 1, Windows 8.x, Windows 10 (professional version); 32 bits. Recommended: Windows 10 Professional 64 bits.
Processor	Minimum: Intel Core i7 at 2.8Ghz.
RAM Memory	Minimum: 4 Gb. Recommended: 8 Gb.
Graphic Card	Minimum: GPU card with minimum 2 Gb of total dedicated VRAM (NVIDIA Card is recommended as it is the most tested).
Hard Disc	5 Gb for the software installation. Additional space may be needed, depending on the volume of cases.
Firewall configuration	Allow access to TCP/IP communication to executables NemoServer.exe and NemoStudio.exe
Network Speed (wired or wireless)	Network local installation. Minimum 1Gb. Internet speed ADSL 50 mb.
Internet Access	Yes, mandatory for licensing control (virtual id).
Mouse	A mouse with three buttons and a mousewheel is highly recommended.
For Apple Macbook, iMac	Bootcamp (native Windows) for the highest performance. Virtual Pc running Windows in a virtual environment. Parallels is recommended. It is advisable to configure the virtual environment with the recommended Windows requirements aforementioned.